



## Developing a high temperature magnetometer to study oceanic basalt magnetisation carriers

Michael Wack, Maxim Alexutin and Jürgen Matzka
Department of Earth and Environmental Sciences, University Munich
Theresienstr. 41, 80333 Munich, MichaelWack@gmx.net

## Why is oceanic basalt such a peculiar magnetisation carrier and what relevance has continuous thermal demagnetisation?




## The "Hotspin 2" - developing an instrument for continuous thermal sensors & Sensitivity demagnetisation of the NRM of rock samples Heating



-35

Sensor arrangement of Hotspin 2. The magnetic sensors are fixed in space, the specimen (cylinder) is moving on a circular path in the y-x-plane.




Magnetic signal (z-component) of sample AGA9\_2 with a magnetization of 55 mA/m measured by Hotspin 2. The noise level is in the order of 10 mA/m.




The Hotspin 2 spinner magnetometer is enclosed by three Helmholtz coils (photo on the left side) to generate a field free environment for the demagnetisation process. Main mechanical features are the rotating sample platform below the circular oven and the sensor holder (sketch on the right side).



References:
[1] U. Bleil, N. Petersen, Variations in magnetization intensity and low-temperature titanomagnetic oxidation of ocean floor basalts, Nature 301 (1983) 384-388.
[2] profiles measured by WOODS HOLE OCEANOGRAPHIC INSTITUTE (CH099L03) and LAMONT-DOHERTY GEOLOGICAL OBSERVATORY (C1604) provided by National Geophysical Data Center (NGDC).
[3] International Geological-Geophysical Atlas of the Atlantic Ocean. Udintsev G.B.(ed). IOC (of UNESCO), Min.Geol.USSR, Ac.Sci. USSR, Moscow, 1989-1990.



The oven encloses the circular sample track and is completely non-magnetic. Hot air is guided through copper pipes which heat the samples by radiation.



Initially, the new hot air heating concept was limited to temperatures up to 140°C (green). After improving the air circulation, 400°C can now be achieved (red).